












SpaceFile®

Interproximal Reduction Kit

Power up your IPR

Versatility and Value

- Swiss-made, diamond-impregnated files
- Uses include breaking contact, reduction, contouring, fine trimming and finishing
- Designed to allow low temperature operation and prevent clogging
- Double-sided or single-sided
- Files, grip and inner tray are autoclavable
- Color coded to match Endodontic file system
- Can be used three ways: Hand Held, With removable grip, in optional air-driven handpiece

DIAMOND-IMPREGNATED		DOUBLE-SIDED	SINGLE-SIDED LEFT	SINGLE-SIDED RIGHT
Assorted Packs		04-500-10	04-500-10L	04-500-10R
Assorted Pack includes (1) Removable Grip, 8 Files (one of each thickness) and 1 Serrated Contact Point Saw				
Single Type Packs				
Polishing A final step in tooth surface restoration.		04-500-60	04-500-60L	04-500-60R
		0.085mm	0.070mm	0.070mm
Finishing For finishing prior to final polish.		04-500-30	04-500-30L	04-500-30R
		0.110mm	0.100mm	0.100mm
XX Fine Pre-polish and finish reduced areas. Also use as initial strip for excessive crowding.		04-500-15	04-500-15L	04-500-15R
		0.140mm	0.123mm	0.123mm
X Fine Pre-polish stripped areas. Also use as initial strip for crowded cases.		04-500-25	04-500-25L	04-500-25R
		0.160mm	0.138mm	0.138mm
Fine Use to initiate reduction in crowded cases. Also use to contour after reduction.		04-500-40	04-500-40L	04-500-40R
		0.270mm	0.180mm	0.180mm
Medium Use for reduction with some crowding or when some space exists. Smaller grain can achieve results before using Coarse strip.		04-500-150	04-500-150L	04-500-150R
		0.330mm	0.210mm	0.210mm
Coarse Primary working strip for IPR. Also use to trim and contour.		04-500-90	04-500-90L	04-500-90R
		0.370mm	0.235mm	0.235mm
X Coarse Use for maximum reduction at start of IPR.		04-500-120	04-500-120L	04-500-120R
		0.490mm	0.290mm	0.290mm
Serrated Contact Point Saw Use to break contact points on enamel, composite, amalgam, porcelain or gold.		04-500-22	—	—
Single-Thickness Packs include Removable Grip and 8 Files				

Removable Grip

When used with strip, provides better hand-held control, resulting in a more precise reduction.



04-500-04
3-Pack



SpaceFile® Contra-Angle Handpiece

- Use with ALL SpaceFile® Files
- Compact design makes reaching posterior regions easy
- Use with most standard slow-speed air motors & compressed air units such as A-dec, KaVo, W&H, NSK, Midwest-E and E-Stylus
- Reciprocating motion of 1.6mm
- Autoclavable
- One-Year Warranty

#40027ANG Contra-Angle Handpiece
#40026MOT Air Motor



Schedule Your Free Lunch and Learn with GAC!

We'll bring lunch to your office and provide a hands-on demonstration for you and your staff!¹

To sign up for your FREE GAC Lunch and Learn, visit:

dentsplysirona.com/orthodontics/ortholearn

1-800-645-5530

¹ Some geographical restrictions may apply.

The Dental Solutions Company

Dentsply Sirona is the world's largest manufacturer of professional dental products and technologies. We create industry-leading positions and platforms across consumables, equipment, technology, and specialty products. We are committed to introducing new, innovative products and complementary, end-to-end solutions to drive better, safer and faster dentistry.

To learn more, visit us online at **dentsplysirona.com**.



Dentsply GAC International
One CA Plaza, Suite 100
Islandia, NY 11749

©2018 Dentsply Sirona
All Rights Reserved
0280-MM-120-500-00 Rev. 03