

# OmniArch® PLUS

## BASE Rx

Dentsply Sirona Orthodontics Base Rx has prescription values that are equivalent to the Roth prescription. ROTH is a registered trademark of Roth Licensing, LLC. All rights reserved.

Maxillary Tooth	T	A	Of	.018		.022	
				Right	Left	Right	Left
U1	12	5	0	J77-111-00	J77-211-00	J77-112-00	J77-212-00
U2	8	9	0	J77-121-00	J77-221-00	J77-122-00	J77-222-00
U3/hk	2	10	0	J77-131-10	J77-231-10	J77-132-10	J77-232-10
U4/univ	-7	0	0	J77-541-00	J77-541-00	J77-542-00	J77-542-00
U4/hk	-7	0	0	J77-141-10	J77-241-10	J77-142-10	J77-242-10
U5/univ	-7	0	0	J77-551-00	J77-551-00	J77-552-00	J77-552-00
U5/hk	-7	0	0	J77-151-10	J77-251-10	J77-152-10	J77-252-10
U4-5/ob	-7	0	0	J77-541-07	J77-541-07	J77-542-07	J77-542-07

Mandibular Tooth	T	A	Of	.018		.022	
				Right	Left	Right	Left
L1-2	-1	0	0	J77-611-00	J77-611-00	J77-612-00	J77-612-00
L3/hk	-11	7	0	J77-331-10	J77-431-10	J77-332-10	J77-432-10
L4/univ	-17	0	0	J77-641-00	J77-641-00	J77-642-00	J77-642-00
L4/hk	-17	0	0	J77-341-10	J77-441-10	J77-342-10	J77-442-10
L4/univ/ob	-17	0	0	J77-641-07	J77-641-07	J77-642-07	J77-642-07
L5/univ	-22	0	0	J77-651-00	J77-651-00	J77-652-00	J77-652-00
L5/hk	-22	0	0	J77-351-10	J77-451-10	J77-352-10	J77-452-10
L5/univ/ob	-22	0	0	J77-651-07	J77-651-07	J77-652-07	J77-652-07

KITS (brackets only)	.018"	.022"
U/L 5-5 Cuspid hk	KJ77-055-13	KJ77-055-23
U/L 5-5 Cuspid Bicuspid hk	KJ77-055-12	KJ77-055-22

## EURO Rx

Euro Rx has prescription values that are equivalent to the \*MBT™ prescription.

Maxillary Tooth	T	A	Of	.018		.022	
				Right	Left	Right	Left
U1	17	4	0	77-111-60	77-211-60	J77-112-60	J77-212-60
U2	10	8	0	77-121-60	77-221-60	J77-122-60	J77-222-60
U3/hk	-7	8	0	77-131-61	77-231-61	J77-132-61	J77-232-61
U4	-7	0	0	77-541-60	77-541-60	J77-542-60	J77-542-60
U4/hk	-7	0	0	77-141-61	77-241-61	J77-142-61	J77-242-61
U5	-7	0	0	77-551-60	77-551-60	J77-552-60	J77-552-60
U5/hk	-7	0	0	77-151-61 7	7-251-61	J77-152-61	J77-252-61

Mandibular Tooth	T	A	Of	.018		.022	
				Right	Left	Right	Left
L1-2	-6	0	0	77-611-60	77-611-60	J77-612-60	J77-612-60
L3/hk	-6	3	0	77-331-61	77-431-61	J77-332-61	J77-432-61
L4	-12	2	0	77-341-60	77-441-60	J77-342-60	J77-442-60
L4/hk	-12	2	0	77-341-61	77-441-61	J77-342-61	J77-442-61
L5	-17	2	0	77-351-60	77-451-60	J77-352-60	J77-452-60
L5/hk	-17	2	0	77-351-61	77-451-61	J77-352-61	J77-452-61

\*MBT is a registered trademark of 3M Unitek. All rights reserved.

KITS (brackets only)	.018"	.022"
U/L 5-5 Cuspid hk		KJ77-055-61
U/L 5-5 Cuspid Bicuspid hk		KJ77-055-62

MicroArch® PLUS  
OmniArch® PLUS  
Next Generation Traditional  
Twin Brackets by GAC



## The Dental Solutions Company

Dentsply Sirona is the world's largest manufacturer of professional dental products and technologies. We create industry-leading positions and platforms across consumables, equipment, technology, and specialty products. We are committed to introducing new, innovative products and complementary, end-to-end solutions to drive better, safer and faster dentistry.

To learn more, visit us online at [dentsplysirona.com](https://dentsplysirona.com).



Dentsply GAC International  
One CA Plaza, Suite 100  
Islandia, NY 11749

©2018 Dentsply Sirona  
All Rights Reserved  
0280-MM-120-J72-03 Rev. 02



MicroArch® PLUS  
OmniArch® PLUS  
Next Generation Traditional  
Twin Brackets by GAC

# Comprehensive Prescription Guide



# It All Adds Up

We weren't looking to make all the other brackets obsolete or impractical. We simply set out to build a better bracket. The PLUS line of products utilize a number of proprietary manufacturing processes that have the ability to deliver increases in torque, angulation and other critical dimensions. When you compare all the advantages we start out with, you can understand why your cases will finish out so well.

Every orthodontic bracket prescription is the sum of its parts. Only PLUS gives you the confidence that comes from an unmatched combination of tangible and intangible benefits. How does your bracket measure up to PLUS?

Dedicated Manufacturing Facility	Yes
Integrated Design Philosophy	Yes
Total Patents Owned	2,500
Annual Clinical Programs	9,572
R&D Team Members	70
Annual New Product Launches	30
Ongoing Clinical Studies	50
Countries Containing Infrastructure	40



## MicroArch® PLUS

### Base Rx

Dentsply Sirona Orthodontics Base Rx has prescription values that are equivalent to the Roth prescription. ROTH is a registered trademark of Roth Licensing, LLC. All rights reserved.

Maxillary				.018		.022	
Tooth	T	A	Of	Right	Left	Right	Left
U1	12	5	0	J72-111-00	J72-211-00	J72-112-00	J72-212-00
U2	8	9	0	J72-121-00	J72-221-00	J72-122-00	J72-222-00
U3/hk	2	10	4	J72-131-10	J72-231-10	J72-132-10	J72-232-10
U4	-7	0	2	J72-141-00	J72-241-00	J72-142-00	J72-242-00
U4/hk	-7	0	2	J72-141-10	J72-241-10	J72-142-10	J72-242-10
U5	-7	0	2	J72-151-00	J72-251-00	J72-152-00	J72-252-00
U5/hk	-7	0	2	J72-151-10	J72-251-10	J72-152-10	J72-252-10
U4-5/ob	-7	0	0	J72-541-08	J72-541-08	J72-542-08	J72-542-08

Mandibular				.018		.022	
Tooth	T	A	Of	Right	Left	Right	Left
L1-2	-1	0	0	J72-611-00	J72-611-00	J72-612-00	J72-612-00
L3/hk	-11	7	2	J72-331-10	J72-431-10	J72-332-10	J72-432-10
L4	-17	0	4	J72-341-00	J72-441-00	J72-342-00	J72-442-00
L4/hk	-17	0	4	J72-341-10	J72-441-10	J72-342-10	J72-442-10
L4/ob	-17	0	0	J72-641-08	J72-641-08	J72-642-08	J72-642-08
L5	-22	0	4	J72-351-00	J72-451-00	J72-352-00	J72-452-00
L5/hk	-22	0	4	J72-351-10	J72-451-10	J72-352-10	J72-452-10
L5/ob	-22	0	0	J72-651-08	J72-651-08	J72-652-08	J72-652-08

KITS (brackets only)		.018"	.022"
U/L 5-5 Cuspid hkm		KJ72-055-13	KJ72-055-23
U/L 5-5 Cuspid Bicuspid hkm		KJ72-055-12	KJ72-055-22
U/L 5-5 Cuspid hk ob			KJ72-055-08

### Euro Rx

Euro Rx has prescription values that are equivalent to the \*MBT™ prescription.

Maxillary				.018		.022	
Tooth	T	A	Of	Right	Left	Right	Left
U1	17	4	0	72-111-60	72-211-60	J72-112-60	J72-212-60
U2	10	8	0	72-121-60	72-221-60	J72-122-60	J72-222-60
U3/hk	-7	8	0	72-131-61	72-231-61	J72-132-61	J72-232-61
U4/univ	-7	0	0	72-541-60	72-541-60	J72-542-60	J72-542-60
U5/univ	-7	0	0	72-551-60	72-551-60	J72-552-60	J72-552-60

Mandibular				.018		.022	
Tooth	T	A	Of	Right	Left	Right	Left
L1-2	-6	0	0	72-611-60	72-611-60	J72-612-60	J72-612-60
L3/hk	-6	3	0	72-331-61	72-431-61	J72-332-61	J72-432-61
L4	-12	2	0	72-341-60	72-441-60	J72-342-60	J72-442-60
L4/hk	-12	2	0	72-341-61	72-441-61	J72-342-61	J72-442-61
L5	-17	2	0	72-351-60	72-451-60	J72-352-60	J72-452-60
L5/hk	-17	2	0	72-351-61	72-451-61	J72-352-61	J72-452-61

\*MBT is a registered trademark of 3M Unitek. All rights reserved.

KITS (brackets only)		.018"	.022"
U/L 5-5 Cuspid hk			KJ72-055-61

## MicroArch® PLUS

### Vertical Slot Rx

The vertical slot, or V-slot, permits use of auxiliaries and biomechanical adjuncts for correction of severely rotated or displaced teeth.


Maxillary				.018		.022	
Tooth	T	A	Of	Right	Left	Right	Left
U1	12	5	0	J72-111-03	J72-211-03	J72-112-03	J72-212-03
U2	8	9	0	J72-121-03	J72-221-03	J72-122-03	J72-222-03
U3/hk	-2	7	0	J72-131-13	J72-231-13	J72-132-13	J72-232-13
U4	-7	0	0	J72-541-03	J72-541-03	J72-542-03	J72-542-03
U4/hk	-7	0	0	J72-141-13	J72-241-13	J72-142-13	J72-242-13
U4-5/ob	-7	0	0	J72-542-33	J72-542-33	J72-542-33	J72-542-33

Mandibular				.018		.022	
Tooth	T	A	Of	Right	Left	Right	Left
L1-2	0	0	0	J72-611-03	J72-611-03	J72-612-03	J72-612-03
L3/hk	-5	7	0	J72-331-13	J72-431-13	J72-332-13	J72-432-13
L4	-11	0	0	J72-641-03	J72-641-03	J72-642-03	J72-642-03
L4/hk	-11	0	0	J72-341-13	J72-441-13	J72-342-13	J72-442-13
L4/ob	-11	0	0	J72-641-33	J72-641-33	J72-642-33	J72-642-33
L4/ob hk	-11	0	0	J72-341-23	J72-441-23	72-342-23	72-442-23
L5	-17	0	0	J72-651-03	J72-651-03	J72-652-03	J72-652-03
L5/hk	-17	0	0	J72-351-13	J72-451-13	J72-352-13	J72-452-13
L5/ob	-17	0	0	J72-651-33	J72-651-33	J72-652-33	J72-652-33

### MicroProgressive Rx

Maxillary				.018	
Tooth	T	A	Of	Right	Left
U1	22	5	0	J72-111-70	J72-211-70
U2	13	9	0	72-121-01	72-221-01
U2/hk	13	9	0	J72-121-11	J72-221-11
U3	7	5	0	72-131-70	72-231-70
U3/hk	7	5	0	J72-131-71	J72-231-71
U4-5	0	0	0	74-541-22	74-541-22
U4-5/hk	0	0	0	74-142-10	74-241-10

Mandibular				.018	
Tooth	T	A	Of	Right	Left
L1-2	-1	0	0	J72-611-00	J72-611-00
L1-2/hk	-1	0	0	72-311-11	72-411-11
L3	7	5	0	72-331-70	72-431-70
L3/hk	7	5	0	J72-331-71	J72-431-71
L4-5	0	0	0	74-641-22	74-641-22
L4-5/hk	0	0	0	74-341-10	74-441-10



Not all skus are available in the PLUS line