



hurrimix™ Automatic Mixer for Alginate and Stone

# Mix with Confidence

The hurrimix™ provides you with highly consistent quality mixes, while reducing mix time by as much as 80% vs hand mixing — a significant savings. The number of rotations achievable with an automatic mixer are not achievable by mixing manually.

1.800.645.5530 [gacintl.com](http://gacintl.com)

# Now You Can Mix Both!

The hurrimix™ increases office efficiencies by standardizing the mixing process to deliver a top quality mix every time.

## Consistent Mix Quality Manual Mixing vs. Automatic Mixing

The hurrimix™ provides you with highly consistent quality mixes, while reducing mix time by as much as 80% vs hand mixing – a significant savings.

The number of rotations achievable with an automatic mixer are not achievable by mixing manually.

Manual mixing can leave bubbles in the mixture and can cause distortion and other imperfections, compromising the quality of the finished work.



Manual Mixing:  
Revolution max 180 rpm



Automatic Mixing with the hurrimix  
Revolution ≈ 600 rpm  
Rotation ≈ 3000 rpm

## Mix Alginate and Stone in 5 Easy Steps



1. Add Powder



2. Add Water



3. Insert Container



4. Start Machine



5. Remove Mix

## Software with 20 Mixing Programs

The hurrimix™ is programmable with 10 options for alginates and 10 programs for stones. Users can modify any programs at any time using its large digital display for simple, rapid operation. You can set programs for specific alginate and stone products resulting in a fast, reliable mixing process, with no time wasted on manually programming the mixer, and guaranteed correct mixing times.

## 2 Year Warranty

The hurrimix™ Automatic Mixer Includes: 1 spatula, 1 pre-mix stick, 1 spatula holder, 4 mixing containers (2 Stone and 2 Alginate), 1 water bottle with a graduated measuring cup, Powercord and DFU.



## Schedule Your FREE Lunch and Learn with GAC

We'll bring lunch to your office and provide a hands-on demonstration for you and your staff.\*

To sign up for your FREE GAC Lunch and Learn, visit:

[dentsplysirona.com/orthodontics/ortholearn](https://dentsplysirona.com/orthodontics/ortholearn)

1-800-645-5530

\*Some geographical restrictions may apply. RTE-019-17  
Dentsply GAC International • One CA Plaza, Suite 100 • Islandia, NY 11749