



imprESSIX®

Color Change & Extra Fast Set Alginate

Because you want to make
a great first impression.


imprESSIX Alginate formulas are tear-resistant, strong and easy to process
helping to ensure that your first impression is a great one.

1.800.645.5530 dentsplysirona.com/orthodontics

Alginate Advantages

imprESSIX Alginates are **stronger & more flexible** than leading competitors*

Tear Strength



*Data on File

Both imprESSIX Formulas:



Easy to remove
from mouth



Dustless



2 casts per
impression



Fills tray
smoothly



Natural
fruit aroma



imprESSIX Color Change Alginate
Tells You When It's Ready.

- Setting time: 60 seconds
- Dimensional stability:
5 days (inside sealed plastic bag)

25-001-09 500g Bag
02-020-09 Case (20) 500g Bags



imprESSIX Extra Fast Set Alginate
Do More in Less Time.

- Setting time: 45 seconds
- Dimensional stability:
2 days (inside sealed plastic bag)

25-001-05 500g Bag
25-020-05 Case (20) 500g Bags

Contact your GAC Rep or call **1.800.645.5530** to order today!

dentsplysirona.com/orthodontics