

CEREC Premium CAM SW 4.5: CEREC MC XL Premium Package, CEREC MC XL, CEREC MC X

Motor set	Left	Right	Combination
Set1 (or Set2**)	Step Bur 12S	Cylinder Pointed Bur 12S	1
Set1 (or Set2**)	Step Bur 12	Cylinder Pointed Bur 12S	2
Set1 (or Set2**)	Step Bur 20	Cylinder Pointed Bur 20	3
Set1	Shaper 25 RZ	Finisher 10	4 * Wet Milling
Set1	Cylinder Bur 12 EF	Cylinder Pointed Bur 12 EF	5 ** Extra Fine
Set1	Shaper 25	Finisher 10	6 *** Dry Milling
Set1 (or Set2**)	Step Bur 20	Cylinder Pointed Bur 12S	7


* Only with grinding/milling units from the following serial numbers: inLab MC XL 129001, CEREC MC XL 129001, CEREC MC XL Premium Package 302001, CEREC MC X 231001 or replacement of the left motor on 2-motor machines or of the left motor in the motor set 1 on 4-motor machines.

** Only for 4-motor grinding/milling units

*** Requires CEREC grinding/milling unit with dry milling configuration and CEREC suction unit.

**** USA/Japan

Description	Figure	Ref.
Step Bur 12S	1	62 40 167
Step Bur 12	2	62 60 025
Cylinder Pointed Bur 12S	3	62 40 159
Cylinder Bur 12 EF	4	65 35 186
Cylinder Pointed Bur 12 EF	5	65 35 178
Step Bur 20	6	62 59 597
Cylinder Pointed Bur 20	7	62 59 589
Shaper 25 RZ	8	64 33 440
Finisher 10	9	62 99 387
Shaper 25	10	62 99 395



Manufacturer	Material	Combination
Sirona	CEREC Blocs C / CEREC Blocs	1 2 5
	CEREC Blocs C PC / CEREC Blocs PC	1 2 5
	CEREC Guide Blocs (for CEREC Guide 2)	4
	CEREC Zirconia	3 4 6
	CEREC Zirconia meso	3 4 6
	inCoris TZI C	3 4 6
	inCoris TZI	3 4 6
	inCoris ZI	3 4 6
	inCoris ZI meso	3 4 6
	inCoris CC	3 4
Dentsply	CELTRA® Duo	1 2 5
3M	Lava™ Ultimate	1 2 5
	Paradigm™ MZ 100 ****	1 2 5
COLTENE	BRILLANT Crios	1 2 5
GC	CERASMART™	1 2 5
	Initial™ LRF	1 5
MERZ DENTAL	artBloc® Temp	3 4
SHOFU	BLOCK HC	1 2 5
	BLOCK HC (2L)	1 2 5

Manufacturer	Material	Combination
VITA	VITABLOCS® MarkII	1 2 5
	VITABLOCS® TriLuxe	1 2 5
	VITABLOCS® TriLuxe forte	1 2 5
	VITA ENAMIC®	1 2 5
	VITA ENAMIC® IS	1
	VITA ENAMIC® multiColor	1 2 5
	VITA SUPRINITY® PC	1 2 5
	VITA SUPRINITY® FC	1 2 5
	VITA CAD-Temp® monoColor	3 4
	VITA CAD-Temp® multiColor	3 4
	VITA CAD-Temp® IS	3 4
	VITA YZ HT	3 4 6
	VITA YZ T	3 4 6
	VITA CAD-Waxx	3 4
	IVOCLAR VIVADENT	IPS Empress® CAD
IPS Empress® CAD Multi		1 2 5
IPS e.max® CAD		1 2 5 7
IPS e.max® CAD (Abutment Solutions)		1
IPS Telio® CAD		3 4
IPS Telio® CAD (Abutment Solutions)		3 4
IPS e.max® ZirCAD LT		3 4 6
IPS e.max® ZirCAD MO		3 4 6

Bestell-Nr. 66 16 150, Ä.-Nr.: 123988