


CEREC SW 4.4.2: CEREC MC

MOTOR SET	LEFT	RIGHT	COMBINATION
Set1	Step Bur 12S	Cylinder Pointed Bur 12S	1
Set1	Step Bur 12	Cylinder Pointed Bur 12S	2
Set1	Step Bur 20	Cylinder Pointed Bur 20	3
Set1	Shaper 25 RZ	Finisher 10	4 * Wet Milling
Set1	Shaper 25	Finisher 10	5 ** Dry Milling

DESCRIPTION	FIGURE	REF.
Step Bur 12S	1	62 40 167
Step Bur 12	2	62 60 025
Cylinder Pointed Bur 12S	3	62 40 159
Step Bur 20	4	62 59 597
Cylinder Pointed Bur 20	5	62 59 589
Shaper 25 RZ	6	64 33 440
Finisher 10	7	62 99 387
Shaper 25	8	62 99 395



* Only with units from the following serial number: CEREC MC 202001

** Requires CEREC milling unit with dry milling configuration and CEREC suction unit.

SIRONA	MATERIAL	COMBINATION
	CEREC Blocs C	1 2
	CEREC Blocs C PC	1 2
	CEREC Blocs C In	1 2
	CEREC Zirconia	3 4 5
	inCoris TZI C	3 4 5
	inCoris TZI	3 4 5

DENTSPLY	MATERIAL	COMBINATION
	CELTRA® DUO	1 2

VITA	MATERIAL	COMBINATION
	VITABLOCS® MarkII	1 2
	VITABLOCS® TriLuxe	1 2
	VITABLOCS® TriLuxe forte	1 2
	VITABLOCS RealLife®	1 2
	VITA ENAMIC®	1 2
	VITA SUPRINITY®	1 2
	VITA SUPRINITY® FC	1 2

IVOCLAR VIVADENT	MATERIAL	COMBINATION
	IPS Empress® CAD	1 2
	IPS Empress® CAD Multi	1 2
	IPS e.max® CAD	1 2

GC	MATERIAL	COMBINATION
	CERASMART™	1 2

3M™ ESPE	MATERIAL	COMBINATION
	Lava™ Ultimate	1 2
	Paradigm™ MZ 100 ***	1 2

SHOFU	MATERIAL	COMBINATION
	BLOCK HC	1 2

COLTENE	MATERIAL	COMBINATION
	BRILLIANT Crios	1 2

*** USA/Japan