

CEREC CAM SW 4.6: CEREC MC XL Premium Package, CEREC MC XL, CEREC MC X

Motoren Set	Links	Rechts	Kombination
Set1 (oder Set2**)	Step Bur 12S	Cylinder Pointed Bur 12S	1
Set1 (oder Set2**)	Step Bur 12	Cylinder Pointed Bur 12S	2
Set1 (oder Set2**)	Step Bur 20	Cylinder Pointed Bur 20	3
Set1	Shaper 25 RZ	Finisher 10	4 ⁽¹⁾ Nassfräsen
Set1	Cylinder Bur 12 EF	Cylinder Pointed Bur 12 EF	5 ⁽²⁾ Extrafein
Set1	Shaper 25	Finisher 10	6 ⁽³⁾ Trockenfräsen
Set1 (oder Set2**)	Step Bur 20	Cylinder Pointed Bur 12S	7

⁽¹⁾ Nur mit Schleif-/ Fräseinheiten ab folgenden Seriennummern: inLab MC XL 129001, CEREC MC XL 129001, CEREC MC XL Premium Package 302001, CEREC MC X 231001 oder Austausch des linken Motors bei 2-motorigen Maschinen oder des linken Motors in Motorset 1 bei 4-motorigen Maschinen.

⁽²⁾ Nur für 4-motorige Schleif-/ Fräseinheiten

⁽³⁾ Setzt CEREC Schleif-/ Fräseinheit mit Trockenfräskonfiguration und CEREC Absaugung voraus.


⁽⁴⁾ USA/Japan

⁽⁵⁾ Weltweit außer Japan

⁽⁶⁾ Japan

⁽⁷⁾ Weltweit außer USA

Bezeichnung	Abbildung	REF.
Step Bur 12S	1	62 40 167
Step Bur 12	2	62 60 025
Cylinder Pointed Bur 12S	3	62 40 159
Cylinder Bur 12 EF	4	65 35 186
Cylinder Pointed Bur 12 EF	5	65 35 178
Step Bur 20	6	62 59 597
Cylinder Pointed Bur 20	7	62 59 589
Shaper 25 RZ	8	64 33 440
Finisher 10	9	62 99 387
Shaper 25	10	62 99 395



Hersteller	Material	Kombination
Sirona	CEREC Blocs C / CEREC Blocs	1 2 5
	CEREC Blocs C PC / CEREC Blocs PC	1 2 5
	CEREC Guide Blocs (für CEREC Guide 2)	4
	CEREC Zirconia	3 4 6
	CEREC Zirconia meso	3 4 6
	inCoris TZI C	3 4 6
	inCoris TZI	3 4 6
	inCoris ZI	3 4 6
	inCoris ZI meso	3 4 6
	inCoris ZI CC	3 4
	Dentsply Sirona	CELTRA® Duo
3M	Lava™ Ultimate	1 2 5
	Paradigm™ MZ 100 ⁽⁴⁾	1 2 5
COLTENE	BRILLIANT Crios	1 2 5
GC	CERASMART™ ⁽⁵⁾	1 2 5
	CERASMART™ 270 ⁽⁶⁾	1 2 5
	CERASMART™ 300 ⁽⁶⁾	1 2 5
	Initial™ LRF	1 5
MERZ DENTAL	artBloc® Temp	3 4
SHOFU	BLOCK HC	1 2 5
	BLOCK HC (2L)	1 2 5
KURARAY NORITAKE	KATANA™ Zirconia Block (ST)	3 4 6
	KATANA™ Zirconia Block (STML)	3 4 6
	KATANA™ Avencia Block ⁽⁶⁾	1 2 5

Hersteller	Material	Kombination
VITA	VITABLOCS® MarkII	1 2 5
	VITABLOCS® TriLuxe	1 2 5
	VITABLOCS® TriLuxe forte	1 2 5
	VITA ENAMIC®	1 2 5
	VITA ENAMIC® IS	1
	VITA ENAMIC® multiColor	1 2 5
	VITA SUPRINITY® PC	1 2 5
	VITA SUPRINITY® FC	1 2 5
	VITA CAD-Temp® monoColor	3 4
	VITA CAD-Temp® multiColor	3 4
	VITA CAD-Temp® IS ⁽⁷⁾	3 4
	VITA YZ HT	3 4 6
	VITA YZ T	3 4 6
	VITA CAD-Waxx	3 4
	IVOCLAR VIVADENT	IPS Empress® CAD
IPS Empress® CAD Multi		1 2 5
IPS e.max® CAD		1 2 5 7
IPS e.max® CAD Abutment		1
IPS e.max® ZirCAD LT		3 4 6
IPS e.max® ZirCAD MO		3 4 6
Telio® CAD		3 4
Telio® CAD Abutment		3 4
Tetric® CAD	1 2 5	

Bestell-Nr. 66 16 143, Ä.-Nr.: 125438